



HOW TO CHANGE DRESSINGS HYGIENICALLY

IN SECONDARY HEALING AND CHRONIC WOUNDS

DRESSING CHANGE IN CHRONIC WOUNDS

PREPARING THE DRESSING CHANGE



Inform the patient about the steps taken in dressing changes, consider pain medication/onset of action



Disinfect work area (disposable gloves!); where applicable, cover work area with sterile barrier impervious to liquids



Hygienic hand disinfection



Readying the required materials



Place ble red wound ensure

Discar

resista

lockin

HANDLING THE WOUND AND DRESSING



Gently remove covering dressing



No-touch technique: Remove wound filler with sterile tweezers or sterile gloves



Selective cleansing of the wound Cleanse aseptic wounds from inside out, cleanse septic wounds from outside in, use one sterile pad/swab per wipe



Directly discard disposable materials/waste in a liner impervious to liquids



Irrigate deep wounds straight from the bottle or with sterile bulb-tip cannula/irrigating catheter



Askina algina straig



with sterile irrigating solution;

follow the manufacturer's

Apply new dressing according to doctor's order

instructions!

Examples from the B. Braun range of products



Prontosan® Wound Gel and Wound Gel X – hydrogels with polihexanide and betaine



Askina® Sorb - calcium alginate with CMC (insert with sterile tweezers)

COMPLETION



Discard disposable and waste materials; disinfect all pointof-care areas



Remove and discard disposable gloves and protective clothing



Perform any other measures ordered, e.g., skin care and application of a compression dressing



Hygienic hand disinfection

Docun

www.woundcare-bbraun.com











patient in most comfortacumbent/sitting position; good lighting conditions



Tie back long hair; wear shortsleeved top and, where neces-I must be easily accessible; sary, wear disposable cap and face mask



Put on disposable apron impervious to liquids; perform hygienic hand disinfection; put on disposable gloves

Change gloves; perform

hygienic hand disinfection



d sharps in puncturent container with g lid



e gloves; perform nic hand disinfection



Deposit reusable instruments right after use in instrument tray



Inspect wound; document wound, possibly



 Calgitrol[®] – calcium
te with silver (apply nt from tube)



nent dressing change



with photograph



Askina® DresSil Border - foam dressing with silicone adhesive base



edges/surroundings with

Askina® Foam - non-adherent foam dressing



Time dressing changes according to condition of wound and the needs and objectives of patient

GOOD TO KNOW

- Do not wear any jewelry (e.g., rings or watches).
- Wear special protective clothing whenever wound is contaminated with resistant organisms.
- Coordinated procedure, train all caregivers involved in patient treatment.
- Always treat septic wounds with the same hygienic diligence as in aseptic wounds.
- Keep the dressing change order: Aseptic wounds first, then contaminated wounds, followed by infected wounds. Care for wounds with resistant pathogens last.
- Dressings in large wounds should be changed by two people (one actually changing, the other assisting), in order to reduce the risk of contaminating the surroundings!
- Never touch wound with bare hands (no-touch technique).
- Ready required materials on work area/dressing trolley:
- non-sterile materials close to patient
- sterile materials away from patient
- Check sterile packaging for expiration date and integrity, and open just before use.
- The healthcare professional uses nonsterile gloves and sterile instruments or sterile gloves on the wound.
- Keep doors and windows closed and do not perform any other activities in the room.



Softa-Man[®] ViscoRub | Composition: 100 ml solution contain: Active ingredients: 45 g Ethanol (100%), 18 g Propanol. Excipients: Purified Water, Butanone, Glycerol, Isopropyl Myristate (Ph.Eur.), Cetearyl Ethylhexanoate, Octyldodecanol (Ph.Eur.), Edetol, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Bisabolol. Therapeutic Indications: Hygienic and surgical hand disinfection. Contraindications: Hypersensitivity (allergy) to Ethanol, Propanol or any of the other ingredients. Side Effects: Cases of local alcohol-induced irritation symptoms (e.g. itching, redness) may occur, especially after frequent application. Moreover, contact allergy is possible. Warnings: Flammable. Keep container tightly closed. Keep away from sources of ignition – No smoking. Avoid contact with eyes. Do not apply on injured skin or mucous membranes. For external use only. Flash point: 21 to 22 °C (DIN 51 755) Last Revision: 02 / 2012 Marketing Authorization Holder: B. Braun Melsungen AG, 34209 Melsungen, Germany

Softasept* N | Composition: 100 g solution contain: Active substances: 74.1 g ethanol and 10.0 g isopropyl alcohol. *Other ingredients*: Purified water. Therapeutic Indications: Skin disinfection before surgical procedures, punctures and injections. **Contraindications**: Hypersensitivity (allergy) to ethanol, isopropyl alcohol. Not suited for antiseptic treatment of mucous membranes or use in the immediate vicinity of the eyes. **Possible Side Effects**: Skin irritation such as redness and burning can occur, especially with frequent use. Contact allergies are also possible. **Warnings**: Highly flammable. Keep container tightly closed. Keep away from sources of ignition – No smoking! Do not spray in open flame! Avoid contact with eyes. Do not use on damaged skin or mucous membranes. For external use only. Flash point 14 *C per DIN 51755. **Marketing Authorisation Holder:** B. Braun Melsungen AG 34209 Melsungen Germany (03/2011)

References

S1 Leitlinie chronische und sekundär heilende Wunden - Hygieneanforderung, erstellt durch: Interdisziplinärer Expertenkonsenz im Arbeitskreis Krankenhaus und Praxishygiene der AWMF (Arbeitsgemeinschaft der wissenschaftlichen medizinischen Fachgesellschaften), 01/2014

Die Anforderungen der Krankenhaushygiene an Wundverband und Verbandwechsel, Deutsche Gesellschaft für Krankenhaushygiene (DGKH) ISSN 2196-5226 (A. Kramer, B. Chergui)

Expertenstandard Pflege von Menschen mit chronischen Wunden, deutsches Netzwerk für Qualitätsentwicklung in der Pflege DNQP

Infektionsprävention in Heimen, Empfehlung der Komission für Krankenhaushygiene und Infektionsprävention beim Robert Koch-Institut (RKI) 2005, Kapitel 6.4.1 Wundverbände

B. Braun Medical AG | Infection Control | Seesatz 17 | CH-6204 Sempach Phone +41 58 258 50 00 | Fax +41 58 258 60 00 | info.bbmch@bbraun.com | www.woundcare-bbraun.com | www.bbraun.com